


Bacteria Filter



King Systems
15011 Herriman Blvd
Noblesville, Indiana, 46060 USA

MPS Medical Product Service
Borngasse 20
35619 Braunsfels, Germany
EC REP

Distributed by:
Respironics Inc.
1001 Murry Ridge Lane
Murrysville, PA 15668 USA



1067702 R15
KGP 10/02/2019

Bacteria Filter

ENGLISH

REF 1014047 Bacteria Filter, Disposable Single/CE

REF 342077 Bacteria Filter Box 10

Intended Use

A bi-directional viral bacteria filter used as an accessory for CPAP and bi-level devices. The filter is single patient disposable, and must be disposed of weekly and between patients.

Caution: U.S. federal law restricts this device to sale by or on the order of a physician.
Note: This filter is not made with natural rubber latex or DEHP.

Compatibility

This accessory is compatible with Philips Respironics BiPAP and CPAP systems.

Note: Some devices may not be available in all markets. For more information, please consult your local Philips Respironics representative.

Warnings:

- Monitor ventilator pressures and replace the blocked filter when pressures no longer meet the desired setting.
- Use duration when used on the inspiratory limb versus the expiratory limb may vary greatly.
- When used on expiratory limb, clogged filter may result in increased expiratory resistance.
- This filter is for single use only and is not intended to be disinfected or cleaned after patient use. Dispose of the filter after patient use. Multi-patient use without disinfection between patients increases the risk of infection.
- This filter is not designed for reuse. Respironics cannot guarantee the performance specifications once this filter has been disinfected or cleaned after patient use.
- Inspect for damage or wear. Discard and replace any components as necessary.

Use Information

The duration of filter use will vary by patient and usage. Monitor for blockage by secretions or other particulate and replace as necessary. The filter usage duration will typically be longer when filtering fresh, dry gas on the inspiratory limb of the circuit than when it is used on the expiratory limb of the circuit due to blockages caused by secretions.

Specifications

Dead Space	66 ml
Bacterial/Viral Filter Efficiency	>99.99%
Resistance (@ 0.5 L/s)	0.7 cm H ₂ O (hPa)/L/s
Male Connector	15 mm I.D./22 mm O.D.
Female Connector	22 mm I.D.

Disposal

Dispose of this device in accordance with local collections and recycling regulations. For more information, visit www.philips.com/recycling.

How to Contact Philips Respironics

Should you experience trouble with this equipment or require assistance, please contact your home care provider. If you need to contact Philips Respironics directly, contact customer service at 1-800-345-6443 (USA or Canada only) or +1-724-387-4000, or go to www.respironics.com to find your local customer service contact information.

Symbols Glossary

Refer to <http://www.symbols.philips.com> for a description of the symbols used on this device and its packaging.

Symbol	Title and Meaning
	Prescription device Caution: U.S. federal law restricts this device by or on the order of a physician.
	Packaging unit To indicate the number of pieces into the package.
	Do not re-use To indicate a medical device that is intended for one use or for use on a single patient during a single procedure.
	Not Made with Natural Rubber Latex Indicates that this product is not made with natural rubber latex.
	Medical Device Indicates that the item is a medical device.
	Unique Device Identifier Indicates the Unique Device Identifier information.
	Date of Manufacture To indicate the date on which a product was manufactured. Country of Manufacturer To indicate the country of manufacture of the product. Note: When applied to the label, "CC" is replaced by the two letter country code.
	Importer Indicates the entity importing the medical device into the locale.